

PRÉIME®

DermaFacial

Préime® DermaFacial
the world's most advanced
aesthetic facial device.



AquaB

Cleanse & Infuse

Spiral Vacuum + Specific Solutions



AquaB

Cleanse and gently exfoliate the skin, while helping to control the size of enlarged pores. This step prepares the skin and infuses nourishing active ingredients that result in lifted and hydrated skin for a fresh and youthful look.

VibroX

Exfoliate & Revive

CO₂ Microdermabrasion + Vibration + Red LED



VibroX

Exfoliate and revive the skin with CO₂ bubbles for a tailored microdermabrasion that penetrates the skin to encourage increased blood circulation. The skin will be cleansed and oxygenated, resulting in reenergized skin appearance.

MicroT

Tighten & Tone

Bipolar Microcurrent 400Hz



MicroT

Tighten and tone the skin to support natural collagen and elastin fiber production. This will help smooth the appearance of the fine lines and wrinkles, while firming the skin and muscle, leaving a youthful look and defined facial contours.

Collagen +

Boost & Rejuvenate

Dual bipolar Radio Frequency 1MHz



Collagen+

Boost and rejuvenate the skin's complexion by inducing controlled thermal penetration, helping stimulate the skin's natural ability to repair itself, stimulating changes in the deep layers of skin allowing the skin's regeneration and creating dense net dermis structure.

UltraB

Nourish & Moisture

Ultrasound 1MHz



UltraB

Nourish and moisturize the skin by assisting transdermal absorption of active ingredients improving moisture retention, absorption of active ingredients contributing to the appearance of healthy skin.